

Transcatheter Mitral Valve Replacement: Current Trends and Future Perspectives

Juan F. Granada, MD

Cardiovascular Research Foundation
New York

Disclosure Statement of Financial Interest

Within the past 12 months, I or my spouse/partner have had a financial interest/arrangement or affiliation with the organization(s) listed below.


Affiliation/Financial Relationship	Company
Institutional Grant/Research Support Co-Founder	Abbott, BSCI, Medtronic, Edwards Cephea Valve Technologies (Abbott)

Natural History of Severe MR in the Real World

1,095 Pts* With 3+/4+ MR and HF (2000 to 2008)

- 74% FMR vs. 21% DMR
- DMR Pts (N=226):
 - 84% SURGERY vs. 16% Med Rx
- FMR Pts (N=814):
 - 36% SURGERY vs. 64% Med Rx
- Un-Operated Med Rx Patients:
 - Lower LVEF
Mean 27% Vs. 42%, $P < 0.0001$
 - Higher STS Score
Median 5.8 Vs. 4.0, $P < 0.001$

Prognosis of Unoperated Patients




* Excluded MVA $\leq 2 \text{ cm}^2$, AR $\geq 2+$, aortic peak velocity $\geq 2.5 \text{ m/s}$, HCM, endocarditis, concomitant AV, Ao or pericardial surgeries, LVAD or OHT.

First Attempt: Transeptal TMVR

CardiAQ
1st Gen.
Steering Snare




CardiAQ
2nd Gen.



Trans-Apical TMVR = Feasibility of Delivery

Multi-Step TMVR = Decrease Catheter Profile




Single-Step TMVR = Procedural Simplification


Early Cephea™ Gen 3 Clinical Experience

South America Feasibility Study


Age	80
Gender	F
MV Disease	4+ SMR
LVEF	32%
LVEDD	5.7 cm
MVA	2.5 cm ²
NYHA	II
STS-PROM	3.3



LOW EF
Hemodynamic Stability



LOW SEPTAL HEIGHT
Depth Control



HIGH PAPILLARY
Precise Delivery


Cephea™ Gen 3: Capsule Separation




Cephea™ Gen 3: Trajectory and Depth Correction




Cephea™ Gen 3: Second Position



Cephea™ Gen 3: Valve Deployment



Cephea™ Gen 3: Final Result



High Screen Failure and Rejection Rates


The Achilles Heel of the TMVR Field

Full D-Shape
Non-Adaptable


D-Shape
Adaptable

Circular
Adaptable

Circular
Non-Adaptable



The Design Concept:
Ideal Anatomical Conditions



The Final Result:
Anatomy-Device Design Mismatch


Why Are TMVR Candidates Rejected?

Phenotypic Clustering Analysis (CHOICE-MI Registry)




2/3 of Screen Failures:
Anatomical Exclusions

Primary MR
Preserved LVEF
MPG >5mmHg
≥moderate MAC
Women
Elevated Age
High TAPSE
High PASP
Prior SAVR




1/3 of Screen Failures:
Clinical Reasons (Frailty)

Secondary MR
High LVESV
High LVEDV
Mixed MR
CAD
Arteriopathy
COPD
≥moderate TR
Hypertension




MR Patients Considered for TMVR


Display Complex Mitral Valve Features




Severe MAC



Multiple Cleft Leaflets



Multiple MR Jets




Bioprosthetic Failures

Special Frame Design Considerations Are Needed!


Severe MAC Prevalence in the SUMMIT Trial

~1/3 of All Screened Patients






Special valve sizes, design considerations, and the use of ancillary tools may be needed to treat these patients

The Risk of LVOT-O Remains High in Specific Patient Subsets (i.e., MAC)



Electrosurgery will play a major role in these patients!

Impact of TA Access in Clinical Outcomes

	TA Tendyne EU Experience ¹ (195 pts)	TA Intrepid ² (252 pts)	TS Intrepid ³ (33 pts)
Procedural Success	94.9%	98.4%	93.9%
Disabling Stroke	4.1%	3.6%	0% 
Major Bleeding	14%	27%	30.9% 
30-Day Mortality	9.3%	13.1%	0% 
1-Year Mortality	28.6%	27.3%	6.7% 
1-Year HFH Rate	25.4%	26.1%	6.7% 
NYHA Class I or II (Baseline vs 1-Year)	22.6% vs 82.5%	19.4% vs 89.1%	30% vs 92%

¹Hell MM. JACC Interv. 2024, ²Bapat V. JACC Interv. 2024, ³Zahr F. JACC Interv. 2023.

EFS TMVR Experience: Procedural Data



	AltaValve (30 pts)	Cephea (24 pts)	Eos (60 pts)	HighLife (30 pts)	Innovalve (45 pts)	Intrepid (33 pts)	M3 (35 pts)
Technical Success	29 (97)	23 (95.8)	60 (100)	27 (90)	43 (95.6)	31 (93.9)	31 (88.6)
In-hospital mortality	0	0	0	0	0	0	0
Device-Related Complications	N/A	1 (4.2)	0	1 (3.3)	0	1 (3.1)	0
LVOT Obstruction	0	2 (8.3)	1 (1.7)	0	0	1 (3.1)	N/A
Conversion to Surgery	1 (3)	1 (4.2)	0	1 (3.3)	N/A	1 (3.1)	0
ASD Closure Performed	N/A	3 (12.5)	2 (3.3)	0	8 (18)	24 (75)	1 (2.9)
Major Access Site Complication	0	0	0	2 (6.7)	0	8 (24.2)	0
≥ Major Bleeding	0	0	N/A	4 (13.3)	3 (6.7)	9 (27.3)	0
Disabling Stroke	0	0	0	0	0	0	N/A

EFS TMVR Experience: 30-Day Results

	AltaValve (30 pts)	Cephea (24 pts)	Eos (60 pts)	HighLife (30 pts)	Innovalve (45 pts)	Intrepid (33 pts)	M3 (35 pts)
All-Cause Mortality	3 (10)	1 (4.2)	3 (5)	3 (10)	2 (4.5)	0	1 (2.9)
CV Mortality	2 (7)	1 (4.2)	3 (5)	1 (3)	N/A	0	N/A
Disabling Stroke	0	0	0	1 (3)	0	0	1/35 (2.9)
≥ MVARC Major Bleeding	N/A	0					0
Major Vascular Complications	N/A	N/A	N/A	3 (10)	0	8 (24.2)	0
NYHA class I or II	29 (96)	13 (65)	41 (68.5)	22 (76)	28 (67)	25 (83)	18 (54.9)
MR ≥2 +	0	0	1 (1.7)	0	1 (2)	0	1 (3.1)
PVL ≥2+		0	N/A	0	N/A	0	N/A
Mean MVG, mmHg	4.2	5.5 ± 2.9	3.8 (2.9-5)	4 (3.3-4.9)	N/A	4.9 (3.7-6.7)	5.6 ± 0.4
MV Reintervention	0	0	1 (1.7)	1 (3)	N/A	1 (3)	2 (5.8)

EFS TMVR Experience: 1-Year Results

	AltaValve* (30 pts)	Eos (60 pts)	HighLife (30 pts)	Intrepid (33 pts)
All-Cause Mortality	4 (14)	12 (20.3)	5 (17)	2 (6.7)
CV Mortality	2 (7)	10 (16.9)	2 (7)	1 (3)
HFH	N/A	16 (27.1)	7 (23)	2 (6.7)
Disabling Stroke	0	1 (1.7)	1 (3)	0
MVARC-Major Bleeding	N/A	18 (30.5)	7 (22)	10 (30.9)
NYHA Class I or II	28 (94)	52 (86.8)	17 (78)	22 (92)
MR ≥ 2	0	0	0	0
PVL ≥ 2	N/A	N/A	1 (4)	0
Mean MVG, mmHg	4.2	3.5 (2.5 - 4.8)	4 (3.3-4.9)	4.6 (3.9 – 5.3)
MV Reintervention	0	3 (5.1)	1 (3)	1 (3)


*6-Month results

TMVR Design Concepts Continue to Evolve

SINGLE STEP (ANNULAR)



INTREPID



CEPHEA




EVOQUE



REVALVE

MULTIPLE STEP (SUB-ANNULAR)




HIGHLIFE



SAPIEN M3



29Fr



INNOVALVE




TIOGA




INNOVHEART

LVOT-O Risk in TMVR is Multi-Factorial





Impact of Frame Design and Anchoring



New Gen. TMVR Frame Designs




Precise Valve Positioning and Axial Correction is Needed in Complex Mitral Cases




Evolution of Adjunctive Interventional Tools Enables the Treatment of Complex Anatomies


Work-Up




A 10.8 cm² – P 125.6 mm
AP 32.6 mm – CC 41.6 mm



Severe MAC





* Surgical Edge-to-Edge Repair




Suture Laceration



TMVR - CEPHEA



Pivotal Trial Timeline* in Major TMVR Programs



 Approved

 Projected Approval

 Pivotal Study

Where Are We in 2025 in the TMVR Field?

- Transfemoral TMVR is a reproducible and safe procedure, with surgical conversion rates below 3%.
- Procedural and long-term outcomes are improving; mortality remains closely linked to patient selection.
- LVOT obstruction remains a key screening barrier; design enhancements and electrosurgery may broaden eligibility
- A single TMVR solution may not suit all; diverse designs are needed to address real-world anatomical variability